

IHE Integration Statement			
Vendor	Product Name	Version	Date
Vital Images, Inc	Vitrea Enterprise Suite	1.2	12/28/2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Access to Radiology Information	Image Display	Multiple Sources	
Scheduled Workflow	Evidence Creator	None	
	Image Display	None	
Internet address for vendor's IHE information: http://www.vitalimages.com/Our_Company/Compliance_Documents.aspx			
Link to Standards Conformance Statements for the Implementation			
DICOM	http://www.vitalimages.com/Our_Company/Compliance_Documents.aspx		
Links to General Information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j	

End of Document